



Giraffe Blue Spot PT Lite Checklist

Customer Information

Customer Name

Service Record #:

System ID or Serial Number

Customer Equipment ID:

Service Task:

Install

Software Version:

Voltage Rating (VAC):

100-240 VAC

Start Date

mm/dd/yyyy

Country

UNITED STATES-US

Modality

MIC

☒ Installation

Visual Inspection

Task	Expected Results	Results	Comments
Power cord check	No signs of damage	<div></div>	<div></div>
Power cord retainer check	Retainer is in place	<div></div>	<div></div>
Overall system parts check	No damaged or missing parts	<div></div>	<div></div>

Task	Expected Results	Results	Comments
Air filter check	No filter damage or dust	<div><div></div><div>▼</div></div>	<div></div>
Labels check	No missing or damaged labels	<div><div></div><div>▼</div></div>	<div></div>
Roll stand check if present	No damaged or missing parts on roll stand	<div><div></div><div>▼</div></div>	<div></div>
Roll stand wheels check if present	wheel locks are functional	<div><div></div><div>▼</div></div>	<div></div>

Installation

Task	Expected Results	Results	Comments
Light Pipe Attachment	Insert the light pipe tip into the light pipe port until the clip is aligned to the top mounting screw hole. Ensure screw mounted properly	<div><div></div><div>▼</div></div>	<div></div>
Mounting the System	Position the system on the accessory dovetail rail. Ensure that dovetail rail tighten by two mounting screws.	<div><div></div><div>▼</div></div>	<div></div>
Check that the light box and dovetail attachment	Ensure light box is securely attached to the accessory dovetail rail or roll stand, i.e. the device does not slide up or down the dovetail rail	<div><div></div><div>▼</div></div>	<div></div>
Verify movement of light pipe	Ensure light pipe moves freely and stays in position	<div><div></div><div>▼</div></div>	<div></div>



Functional Checks

Task	Expected Results	Results	Comments
Lighpipe positioning check	Lighpipe moves freely and stays in position	<div><div></div><div>▼</div></div>	<div></div>
Turn Hour meter and check its funcion	Confirm that the hour meter display turns on and the standby/on indicator illuminates in green	<div><div></div><div>▼</div></div>	<div></div>
Standby/on switch check	Ensure Press and release the standby/on switch, and confirm the over-temperature/under-temperature red indicator flashes briefly and the standby/on indicator changes to blue.	<div><div></div><div>▼</div></div>	<div></div>
Fan check	Fan is operating	<div><div></div><div>▼</div></div>	<div></div>

Task	Expected Results	Results	Comments
Light check	Therapeutic light is on.	<div><div></div></div>	<div></div>

Light Intensity Check

Task	Expected Results	Observed Results	Results	Comments
Light Focus diameter	Ensure light Focus spot diameter is 35.5 cm (14") or when the shade is 38 cm (15") from the bed surface. Ensure reading measured after 5 min of bulb switched ON	-	<div><div></div></div>	<div></div>
5-point light measurment	Average spectral irradiance is at least 27 $\mu\text{W} \cdot \text{cm}^{-2} \cdot \text{nm}^{-1}$	<div>Spectral Irradiance Measurements: Point 1 (Center):($\mu\text{W} \cdot \text{cm}^{-2} \cdot \text{nm}^{-1}$) <div></div> Point 2 ($\mu\text{W} \cdot \text{cm}^{-2} \cdot \text{nm}^{-1}$) <div></div> Point 3 ($\mu\text{W} \cdot \text{cm}^{-2} \cdot \text{nm}^{-1}$) <div></div> Point 4 ($\mu\text{W} \cdot \text{cm}^{-2} \cdot \text{nm}^{-1}$) <div></div> Point 5 ($\mu\text{W} \cdot \text{cm}^{-2} \cdot \text{nm}^{-1}$) <div></div> Average: ($\mu\text{W} \cdot \text{cm}^{-2} \cdot \text{nm}^{-1}$) <div></div></div>	<div><div></div></div>	<div></div>

Task	Expected Results	Observed Results	Results	Comments
Lowest to Highest measurement ratio	The lowest to highest measurement ratio shall be more than 0.4	Lowest to Highest Ratio: <div></div>	<div>▼</div>	<div></div>

Electrical Safety Test

Task	Expected Results	Observed Results	Results	Comments
Ground resistance check	<div>◇ Unit plugged to wall outlet: measure the resistance between power cord ground terminal and the potential equalization stud on the systemGround resistance is less than 0.1 ohm</div> <div>◇ Unit plugged to accessory outlet: measure the resistance between the power cord ground terminal of the microenvironment and the potential equalization stud on the Giraffe Blue Spot PT Lite. Ground resistance is less than 0.2 ohm</div>	<div>◇ Unit plugged to wall outlet: Ground resistance: (ohm) <div></div></div> <div>◇ Unit plugged to accessory outlet: Ground resistance: (ohm) <div></div></div>	<div>▼</div>	<div></div>

PT device plugged to: ☐ Wall Outlet
☐ Accessory Outlet

Comments, tools, and attachments

Comments Character(s) 2000

Tools used

Description	Serial Number	Bar Code	Cal Due Date	
<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text" value="mm/dd,"/>	<input type="button" value="ADD"/>

Attachments

No file chosen

GE Representative

FSEs completing this PM

▼

FSE Name	SSO #	SR #
<input type="text"/>	<input type="text"/>	<input type="text"/>

Date complete

Customer Signature

☐ Yes ☒ Not Required

SUBMIT FORM

SAVE CURRENT FORM

VALIDATE FORM

CLEAR FORM

Document Number: BBCZH1165 Revision: 3